

# PRODUCT CATALOGUE

## Maxillo Facial Implants



**UNIGENERATION SDN BHD**

21-2 Jalan Putra Mahkota 7/7B, Putra Heights, 47650 Subang Jaya, Selangor, Malaysia

Tel : +603 5614 1599 Fax : +603 5614 5339 Email : [contact@unigeneration.com.my](mailto:contact@unigeneration.com.my)

[www.unigeneration.com.my](http://www.unigeneration.com.my)



Orthognathic titanium plate:  
 1) 0.8mm thickness;  
 2) 4holes;  
 3) L-shape (bevel);  
 4) long stem



Cranio-maxillofacial Osteosynthesis system:  
 Mandible titanium plate:  
 1) 1.0mm thickness;  
 2) 23.6mm length;  
 3) 16.6mm width;  
 4) 5 holes;  
 5) L-shape (bevel);  
 6) long bridge.




Mandible titanium plate:  
 1) 1.0mm thickness;  
 2) 23.6mm length;  
 3) 10.6mm width;  
 4) 4 holes;  
 5) L-shape (bevel);  
 6) Long stem



Mandible titanium plate:  
 1) 1.0mm thickness;  
 2) 22.1mm length;  
 3) 14.6mm width;  
 4) 4 holes;  
 5) T-shape;




Mandible titanium plate:  
 1) 1.0mm thickness;  
 2) 24.7mm length;  
 3) 13.1mm width;  
 4) 6 holes;  
 5) X-shape;



Mandible titanium plate:  
 1) 1.0mm thickness;  
 2) 22.6mm length;  
 3) 12.6mm width;  
 4) 4 holes;  
 5) Y-shape;  
 5) short stem.




Mandible titanium plate:  
 1) 1.0mm thickness;  
 2) 28.6mm length;  
 3) 12.6mm width;  
 4) 4 holes;  
 5) Y-shape;  
 6) long stem;



Mandible titanium plate:  
 1) 1.0mm thickness;  
 2) 4holes;  
 3) Straight;



Mandible titanium plate:  
 1) 1.0mm thickness;  
 2) 2 holes;  
 3) Straight;  
 4) Stem.



1) 1.0mm thickness;  
 2) 6holes;  
 3) straight;  
 4) short stem;



0.6mm thickness;  
 11.8mm length;  
 7.8mm width;  
 4 holes;  
 L-shape, short




Maxilla fixation plate Specification:  
 0.6mm thickness;  
 16.3mm length;  
 7.8mm width;  
 4 holes;  
 L-shape  
 Short stem;




0.6mm thickness;  
16.3mm length;  
3.8mm width;  
5 holes;  
T-shape.



0.6mm thickness;  
18.7mm length;  
3.8mm width;  
4 holes;  
Stem;  
V-shape (arc);  
Short bridge



0.6mm thickness;  
22.8mm length;  
3.8mm width;  
4 holes;  
Stem;  
V-shape (arc);  
Long bridge.



0.6mm thickness;  
26.3mm length;  
3.8mm width;  
5 holes;  
V-shape (arc);  
Stem



0.6mm thickness;  
26.3mm length;  
3.8mm width;  
5 holes;  
V-shape (arc);  
Stem



0.6mm thickness;  
30.1mm length;  
3.8mm width;  
8 holes;  
V-shape (arc)



0.6mm thickness;  
36.1mm length;  
3.8mm width;  
10 holes;  
V-shape (arc);



0.6mm thickness;  
18.1mm length;  
9.3mm width;  
6 holes;  
X-shape



0.6mm thickness;  
18.1mm length;  
12.8mm width;  
4 holes;  
X-shape;



0.6mm thickness;  
18.1mm length;  
12.8mm width;  
5 holes;  
X-shape.




0.6mm thickness;  
20.3mm length;  
12.8mm width;  
4 holes;  
Y-shape;  
Big



1) 0.6mm thickness  
2) 6holes  
3) straight, stem




- 1) 0.8mm thickness;
- 2) 4holes;
- 3) L-shape (right);
- 4) middle stem




- Orthognathic titanium plate:
- 1) 0.8mm thickness;
  - 2) 4holes;
  - 3) L-shape (bevel);
  - 4) long stem




- Orthognathic titanium plate:
- 1) 0.8mm thickness;
  - 2) 4holes;
  - 3) L-shape (bevel);
  - 4) long stem




- Orthognathic titanium plate:
- 1) 0.8mm thickness;
  - 2) 6holes;
  - 3) Straight;




- Orthognathic titanium plate:
- 1) 0.8mm thickness;
  - 2) 6holes;
  - 3) Straight;




- Orthognathic titanium plate:
- 1) 0.8mm thickness;
  - 2) 8 holes;
  - 3) Straight;



Orthognathic titanium plate:  
 1) 0.8mm thickness;  
 2) 4holes;  
 3) straight;  
 4) stem



Titanium mesh plate:  
 1) 0.6mm thickness;  
 2) Size: 100mmx100mm  
 3) Two dimension



Titanium plate for genioplasty:

YX0.8-02.1:  
 0.8mm thickness;  
 16.4mm length;  
 4.8mm width;  
 4 holes;  
 X-shape;  
 4mm spacing;

YX0.8-02.2:  
 0.8mm thickness;  
 16.4mm length;  
 4.8mm width;  
 4 holes;  
 X-shape;  
 6mm spacing;

YX0.8-02.3:  
 0.8mm thickness;  
 16.4mm length;  
 4.8mm width;  
 4 holes;  
 X-shape;  
 8mm spacing;

YX0.8-02.4:  
 0.8mm thickness;  
 16.4mm length;  
 4.8mm width;  
 4 holes;  
 X-shape;  
 10mm spacing