

PRODUCT CATALOGUE

ARTHROSCOPY




UNIGENERATION SDN BHD

21-2 Jalan Putra Mahkota 7/7B, Putra Heights, 47650 Subang Jaya, Selangor, Malaysia

Tel : +603 5614 1599 Fax : +603 5614 5339 Email : contact@unigeneration.com.my

www.unigeneration.com.my



Loop Endobutton


The Loop endobutton product family was designed for femoral fixation in cruciate ligament reconstructions for hamstring and BTB grafts. Each Loop endobutton ensures a reproducible and reliable fixation, regardless of the quality of the cancellous bone. The tactile feedback flipping allows safe placement of LOOP ENDOBUTTONS. An additional lateral incision is not necessary.

- Consistent cortical fixation
- Tactile flipping Feedback


Interference Screw




Fully Thread Interference Screw



Cross Pin







Tibial Post Fixation Screw



Ancor Preloaded



Tibia Base Plate



Stapel

